

PCT

WORLD INTELLECTUAL PROPERTY ORGANIZATION  
International Bureau



INTERNATIONAL APPLICATION PUBLISHED UNDER THE PATENT COOPERATION TREATY (PCT)

<b>(51) International Patent Classification <sup>5</sup> :</b> <b>G01N 33/72, 33/96</b>	<b>A1</b>	<b>(11) International Publication Number:</b> <b>WO 90/11527</b> <b>(43) International Publication Date:</b> 4 October 1990 (04.10.90)
<b>(21) International Application Number:</b> PCT/US90/01422 <b>(22) International Filing Date:</b> 16 March 1990 (16.03.90)  <b>(30) Priority data:</b> 328,622 27 March 1989 (27.03.89) US  <b>(60) Parent Application or Grant</b> <b>(63) Related by Continuation</b> US 328,622 (CIP) Filed on 27 March 1989 (27.03.89)  <b>(71) Applicant (for all designated States except US):</b> BIONOS- TICS, INC. [US/US]; Two Craig Road, Acton, MA 01720 (US).	<b>(72) Inventor; and</b> <b>(75) Inventor/Applicant (for US only) :</b> CHIANG, Ching [US/ US]; 8 Fernwood Road, Acton, MA 01720 (US).  <b>(74) Agents:</b> BROOK, David, E. et al.; Hamilton, Brook, Smith & Reynolds, Two Militia Drive, Lexington, MA 02173 (US).  <b>(81) Designated States:</b> AT (European patent), BE (European + patent), CA, CH (European patent), DE (European pa- tent), DK (European patent), ES (European patent), FR (European patent), GB (European patent), IT (European patent), JP, LU (European patent), NL (European pa- tent), SE (European patent), SU, US.	
<b>Published</b> <i>With international search report.</i>		
<b>(54) Title:</b> A STABLE HEMOGLOBIN REFERENCE SOLUTION  <b>(57) Abstract</b>  A stable reference solution for calibrating and monitoring blood gas instrumentation is disclosed. The solution comprises an aqueous mixture containing a hemoglobin solution derived from a mammalian source which comprises at least about 95 % reduced hemoglobin. The solution additionally contains a bicarbonate buffer and a metal catalyzed methemoglobin reducing system and an organic buffer. To provide a control element having a variety of properties similar to fresh, whole, human blood, the reference solution is stored in a sealed ampule under an inert atmosphere containing CO <sub>2</sub> until just prior to use. The ampule is subsequently opened and the solution is equilibrated with a gas mixture having components which provide gasses in a physiological range to the solution. The solution is characterized by exhibiting properties similar to fresh blood subsequent to equilibration and an extended storage life prior to equilibration.		

### DESIGNATIONS OF "DE"

Until further notice, any designation of "DE" in any international application whose international filing date is prior to October 3, 1990, shall have effect in the territory of the Federal Republic of Germany with the exception of the territory of the former German Democratic Republic.

#### *FOR THE PURPOSES OF INFORMATION ONLY*

Codes used to identify States party to the PCT on the front pages of pamphlets publishing international applications under the PCT.

AT	Austria	ES	Spain	MG	Madagascar
AU	Australia	FI	Finland	ML	Mali
BB	Barbados	FR	France	MR	Mauritania
BE	Belgium	GA	Gabon	MW	Malawi
BF	Burkina Faso	GB	United Kingdom	NL	Netherlands
BG	Bulgaria	HU	Hungary	NO	Norway
BJ	Benin	IT	Italy	RO	Romania
BR	Brazil	JP	Japan	SD	Sudan
CA	Canada	KP	Democratic People's Republic of Korea	SE	Sweden
CF	Central African Republic	KR	Republic of Korea	SN	Senegal
CG	Congo	LI	Liechtenstein	SU	Soviet Union
CH	Switzerland	LK	Sri Lanka	TD	Chad
CM	Cameroon	LU	Luxembourg	TG	Togo
DE	Germany, Federal Republic of	MC	Monaco	US	United States of America
DK	Denmark				

-1-

A Stable Hemoglobin Reference SolutionBackground of the Invention

In practice, the use of analytical equipment for measuring various parameters of blood requires utilization of control solutions which exhibit properties as close as possible to those observed in fresh normal human blood. By measuring known parameters of control solutions, the instrumentation can be monitored and calibrated to allow highly

10 accurate measurements of patient blood samples.

One approach to monitoring the reliability and accuracy of instruments that measure partial pressures of  $\text{CO}_2$  and  $\text{O}_2$  in blood is with the use of samples of fresh human blood which has been tonometered with gas mixtures having known amounts of  $\text{CO}_2$ ,  $\text{O}_2$  and  $\text{N}_2$ . Such a process is described by Burnett in Clinical Chemistry, 27(10):1761 (1981).

When the blood has been properly tonometered, the sample will have precise and fixed partial pressures of  $\text{CO}_2$  and  $\text{O}_2$ . These prepared samples can then be introduced into the analytical instrument and  $\text{pCO}_2$  and  $\text{pO}_2$  values determined. Since the sample is similar to a patient blood specimen, but with known  $\text{pCO}_2$  and  $\text{pO}_2$  values, the instrument can be considered reliable for measuring unknown patient samples if the test values of the tonometered blood equal the theoretical values based on the gas mixture used for tonometry.

- 2 -

Although the use of tonometered blood is considered to be satisfactory for monitoring blood gas instrumentation, the approach has a number of drawbacks which limit its use in all but a very  
05 small percentage of laboratories.

For example, blood samples derived from human sources are susceptible to infectious agents, including hepatitis virus and HIV which can pose serious health hazards to laboratory personnel who  
10 must perform the tonometry and testing of the sample. In addition, the instruments commonly used to measure  $p\text{CO}_2$  and  $p\text{O}_2$  also measure blood pH. Since the tonometered blood does not have a known pH value, the sample cannot be used for monitoring the  
15 pH measurements and a separate pH control standard must be used.

Similarly, many laboratories that perform measurements of  $p\text{CO}_2$ ,  $p\text{O}_2$  and pH on blood samples also measure total hemoglobin and hemoglobin fractions on a CO-Oximeter that is located near the  
20 blood gas instrument. However, since the blood used for tonometry is acquired from random patient samples, the tonometered sample has no known hemoglobin value, and therefore it is not useful for  
25 monitoring the CO-Oximeter. Consequently, a separate control standard is required for this instrument also.

Finally, the entire procedure for properly preparing tonometered blood samples requires dis-  
30 ciplined techniques and many laboratories lack

-3-

trained personnel as well as the time to prepare the samples.

Because of these disadvantages, most laboratories use control standards which mimic human blood but have properties quite different than fresh blood. For example, buffered aqueous solutions which have been tonometered with  $\text{CO}_2$  and  $\text{O}_2$  are often used. These materials are assayed for predetermined values for pH,  $\text{pCO}_2$  and  $\text{pO}_2$ . However, in composition, physical properties and chemical properties they differ greatly from whole blood.

Other control standards comprise buffered suspensions of modified human red blood cells or hemoglobin solutions prepared from lysed red blood cells. These materials have some properties which more closely approximate actual blood than do the aqueous based controls, but their  $\text{pO}_2$  buffering action and inability to provide the  $\text{O}_2$  saturation characteristics of fresh blood, cause these materials to perform more like aqueous solutions than tonometered fresh blood. Furthermore, since they are prepared from human blood, the health risk to technicians is not eliminated.

In summary, the commercial blood gas controls which are used instead of tonometered fresh human blood are generally considered a compromise between convenience, economy and the ideal control standard.

A need exists for a reference solution which is not susceptible to infection, which can be packaged to eliminate the sample collection and preparation steps necessary when blood samples are used and is

-4-

stable in the packaged form. Further, the solution should be capable of being tonometered in the same manner as fresh human blood to provide a control standard that has O<sub>2</sub> saturation characteristics and other properties similar to fresh human blood, but unlike the blood samples, will have predetermined pH and hemoglobin values for monitoring the instrument performance in the measurement of these parameters as well as pO<sub>2</sub> and pCO<sub>2</sub>.

Although some of the commercial blood based materials can be used as prepackaged fluids for tonometry, the inability of these materials to be manufactured and stored without the oxidation of a significant percent of the hemoglobin to methemoglobin (which does not bind with oxygen), causes the solution to lose the oxygen saturation properties of fresh blood, and therefore makes the material unsuitable as a substitute for blood as a tonometry solution.

Attempts to provide a system for reducing methemoglobin content in a blood-based material have been the subject of a variety of scientific studies. For example, in U.S. Patent 4,485,174 to Chiang et al., a "methemoglobin reductase" enzyme system is described as a means for maintaining a low methemoglobin level. This system, however, has demonstrated only limited usefulness, since the supply of methemoglobin reducing reagents can be exhausted when the hemoglobin solution is stored under an oxygen containing atmosphere.

- 5 -

Other blood-gas control solutions and methods have been described in U.S. Patent 3,859,049; U.S. Patent 3,973,913; U.S. Patent 4,001,142; and U.S. Patent 4,469,792. Additionally, a blood-gas control  
05 solution and method has been described by Steiner et al., Clinical Chemistry, 24, 793 (1978). Each of these, however, described a control standard having a limited storage life and/or a chemical formulation which provides physiologically inaccurate values.

10 Thus, a need still exists for a blood-based reference solution which, after equilibration with an appropriate gas mixture, can be used to monitor blood-gas analysis equipment. The reference solution should preferably have an extended storage  
15 life, provide physiologically accurate blood-gas values, and provide uniform values among a large number of samples.

#### Summary of the Invention

The present invention pertains to a stable  
20 hemoglobin solution which can be used as a control reference for blood-gas analysis instrumentation as well as CO-Oximeters. More specifically, the present invention pertains to an aqueous hemoglobin-containing solution which can be equilibrated with a  
25 variety of gasses generally present in blood to provide a solution having a variety of parameters similar to those parameters normally found in fresh whole blood. Such an equilibrated solution has utility as a control in blood-gas analysis and  
30 CO-Oximetry instrumentation.

-6-

One embodiment of the invention is an aqueous solution comprising:

- a) hemoglobin;
- 05 b) a source of bicarbonate ions to stabilize the pH of the solution following equilibration with a predetermined gas mixture;
- c) a source of thio functional groups to maintain a low level of methemoglobin;
- 10 d) a source of metal ions to catalyze the methemoglobin reduction reaction;
- e) at least one organic polyphosphate; and
- f) an organic buffer.

To maintain anaerobic conditions, the mixture is stored under an inert atmosphere (preferably  
15 nitrogen) containing CO<sub>2</sub>. The organic polyphosphates act to stabilize the hemoglobin molecules in a spatial configuration which is not easily amenable to oxidation to methemoglobin.

Additionally, the solution can contain a  
20 variety of additional components, including, but not limited to, antibacterial and antifungal agents and red blood cell lysing agents.

The solution is maintained in a sealed environment under anaerobic conditions until shortly before  
25 use. When its use as a control is desired, the solution can be equilibrated with physiological gasses via a method such as tonometry. The level of methemoglobin remains low due to the presence of the thio groups and metal ion catalyst which act to  
30 catalytically reduce any methemoglobin which might



- 7 -

form during tonometry. Once equilibrated, the system provides values of pH,  $pCO_2$  and  $pO_2$  suitable for monitoring Blood Gas/pH analysis instrumentation. The solution also provides physiological levels of hemoglobin and fractions of hemoglobin species for calibrating and monitoring CO-Oximetry instrumentation. The solution can also be used as a reference solution for electrolytes ( $Na^+$ ,  $K^+$ ,  $Ca^{++}$ ,  $Cl^-$ ). Finally the solution also exhibits oxygen saturation characteristics similar to those of fresh whole blood.

#### Detailed Description of the Invention

The present invention is a stable reference solution for, inter alia, blood gas and CO-Oximetry instrumentation. Unlike solutions previously employed for calibration and monitoring of blood gas/CO-Oximetry equipment in which physiological gasses are present in solution during storage, the present solution comprises a stable, aqueous, hemoglobin solution which is equilibrated with physiological gasses just prior to use. This solution, and equilibration method associated therewith, provide a blood gas reference solution having a storage life which is longer than that of previously employed controls. This is achieved in part through the use of a metal ion catalyzed methemoglobin reducing system which serves to inhibit formation of undesirable methemoglobin during storage, gas equilibration and use.

- 8 -

A cell-free hemoglobin solution is prepared from healthy animal blood. In the preferred embodiment, the hemoglobin source is bovine, rather than human, blood. Although modified human red blood cell solutions, as well as stroma-free human hemoglobin solutions, are commercially available, these human-source solutions are susceptible to contamination by a variety of human infectious agents including the hepatitis virus and HIV. The use of non-human blood products in blood-based reference solutions is therefore desirable. As such, the use of bovine blood collected from healthy animals is preferred due to its similarity to human blood in both physical and chemical characteristics.

Additionally, the present invention exhibits oxygen saturation characteristics which are more similar to those of fresh whole blood than have been seen in other control standards. For example, when oxygen saturation is plotted against partial pressure, known reference solutions exhibit a hyperbolic curve. In contrast, the present solution, like fresh whole blood, exhibits a sigmoid curve. Furthermore, a value for the required oxygen partial pressure to achieve 50% saturation, the so-called  $P_{50}$  value, has been determined for the reference solution of this invention. This  $P_{50}$  values, obtained for tonometered solutions made from bovine erythrocyte derivatives, were in the range of 18 to 34 mmHg, with a nominal value of 26. This value is very close to that obtained from fresh human blood, generally about 27 mmHg.

- 9 -

The solution of red blood cells is repeatedly washed with saline. Subsequently, the washed red blood cells are combined with deionized water. At this point, some red blood cell lysis may occur.

05 Next, the remaining red blood cells are lysed by the addition of a lysing agent. In the preferred embodiment, the lysing agent is benzethonium chloride. This compound is selected because it also provides antibacterial and antifungal activity to

10 the solution in addition to acting as a cell lysing agent. Optionally at this point, an additional compound, such as the preferred 2-phenoxyethanol, can be added to the solution to provide additional antibacterial and antifungal activity to the

15 solution. The resulting solution contains completely lysed red blood stroma and hemoglobin in solution. The stroma can be removed at this point by centrifugation or it can be removed after subsequent processing steps.

20 The methemoglobin reducing system is now added to the completely lysed hemoglobin solution. The preferred methemoglobin reducing system comprises a sulfhydro compound selected from the class of chemicals which have an effective thio (-SH) functional group. Such chemicals include, but are not

25 limited to the reduced form of glutathione, mercaptoethanol, cysteine and their derivatives. A catalytic amount of a metal ion is provided to the solution via salts of selenium, tellurium, and

30 copper which are added to the solution to catalyze

-10-

the methemoglobin reduction reaction. Such a metal ion catalyzed methemoglobin reducing system is described by Masukawa and Iwata in Life Sciences, 21(5), 695 (1977), the teachings of which are  
05 incorporated herein by reference. In this reference, selenite, selenate, and selenocystine are described as catalysts for reactions in which methemoglobin is reduced by glutathione, 2-mercaptoethylamine and cysteine.

10 ~~Additionally, further resistance to methemo-~~  
globin formation can be imparted to the solution by the addition of a polyphosphate compound. This compound stabilizes the hemoglobin molecules in a spatial configuration which is not readily oxidized  
15 to the methemoglobin form. Preferred polyphosphate compounds include 2,3-diphosphoglycerate, inositol hexaphosphate (also referred to as phytic acid), and adenosine triphosphate. The use of a polyphosphate to stabilize hemoglobin was described in greater  
20 detail by Mansouri in Hemoglobin, 5(6), 579 (1981), the teachings of which are incorporated herein by reference.

Following the addition of the hemoglobin stabilizing systems described above, the pH of the  
25 solution is adjusted to within a physiological range. In the preferred embodiment, the pH is adjusted to between about 7.0 and 7.6 using aqueous NaOH under an inert atmosphere such as nitrogen or noble gas. A source of bicarbonate ion, preferably  
30 sodium bicarbonate and a suitable organic buffer selected from Good's buffers and their derivatives

-11-

are added to buffer the pH. Preferable buffers include N-2-hydroxyethylpiperazine-N'-2-ethane sulfonic acid (HEPES;  $pK_a = 7.31$ ), 3-[N-bis(hydroxyethyl)-amino]-2-hydroxypropane sulfonic acid (DIPSO; 05  $pK_a = 7.35$ ), piperazine-N-N'-Bis (2-hydroxypropane sulfonic acid) (POPSO;  $pK_a = 7.63$ ), 3-[N-(Tris-hydroxymethyl)methylamino]-2-hydroxypropane sulfonic acid (TAPSO;  $pK_a = 7.39$ ), and N-hydroxyethyl-10 piperazine-N'-2-hydroxypropane sulfonic acid (HEPPSO;  $pK_a = 7.73$ ).

The solution is passed through a series of membrane filters having a final filter of a maximum pore size of 0.2  $\mu$ m to remove any red blood cell stroma which remains. Optionally, an antifoaming 15 agent can be added at this point.

The hemoglobin solution is then further bubbled with nitrogen or other inert gasses which contain between about 4 to about 7% carbon dioxide. This nitrogen bubbling step serves to convert any oxy- 20 hemoglobin in solution to reduced hemoglobin and to maintain anaerobic conditions. The carbon dioxide in the gas mixture is used to maintain an acceptable pH value between about 7.0 and 7.2 in the solution. The solution is then sealed in a container, such as 25 glass ampule, under an atmosphere like that bubbled through the solution.

The resulting hemoglobin contained in the aqueous solution preferably comprises hemoglobin fractions of at least about 95% reduced hemoglobin, 30 less than 3% methemoglobin, less than 2% oxyhemoglobin, less than 2% carboxyhemoglobin and less than

-12-

0.3% by volume oxygen. The solution has a purple-blue color as a result of the reduced hemoglobin contained therein. As there are few, if any, oxidizing agents (including molecular oxygen) in solution, the formation of methemoglobin is almost completely eliminated. As there is only negligible methemoglobin in solution, the methemoglobin reducing system is not utilized at this point. Rather, the system will serve to prevent methemoglobin

production during a subsequent gas equilibration just prior to solution use.

When the solution is to be used for monitoring of blood gas/CO-Oximetry instruments, it must first be equilibrated with a physiological gas mixture. In the preferred method, this equilibration is performed via tonometry. Tonometry is the term used to describe the process whereby a liquid is brought into equilibrium with a gas phase by bubbling the desired gas through the liquid or by forming a liquid film with large surface area interacting with the incoming gas. A detailed discussion of tonometry is provided by R.W. Burnett in Clinical Chemistry, 27(10), 1761 (1981), the teachings of which are incorporated herein by reference. During the tonometry (or even when exposed to ambient air) the solution rapidly changes from the purple-blue color indicative of the presence of reduced hemoglobin to a red color indicative of the presence of oxyhemoglobin. The amount of methemoglobin produced during the tonometry remains minimal due to the metal ion catalyzed methemoglobin reduction system described previously.

-13-

Upon completion of the equilibration, the solution will contain known, predetermined levels of pH,  $p\text{CO}_2$  and  $p\text{O}_2$ . Thus the solution can be analyzed by a Blood Gas/pH Analyzer and can be used to

05 monitor the instrumentation and ensure a minimum drift from the set-point. The solution can also be used for the calibration and monitoring of CO-Oximetry equipment used to measure levels of total hemoglobin and fractions of hemoglobin

10 species.

In a preferred embodiment, the solution is representative of normal human blood. In this embodiment, the solution has the following composition:

15	<u>Component</u>	<u>Range</u>	<u>Specific</u>
	Total Bovine		
	Hemoglobin	13-15 g/dl	14 g/dl
	2-Phenoxyethanol	0.09-0.11 g/dl	0.1 g/dl
	Benzethonium Chloride	0.09-0.11 g/dl	0.1 g/dl
20	Sodium Selenite	1-20 uM	2 uM
	Phytic Acid	0.1-1.4 mM	0.14 mM
	Reduced Glutathione	0.5-8 mM	0.8 mM
	Sodium Bicarbonate	28.5-30.5 mM	29.16 mM
	HEPES	20-40 mM	25 mM.

25 Upon use, this solution is equilibrated via tonometry with a gas mixture which comprises by volume 5%  $\text{CO}_2$ , 12%  $\text{O}_2$  and 83%  $\text{N}_2$ . This equilibration is carried out at 37°C until the pH is within the range of about 7.38 to about 7.42, the  $p\text{CO}_2$  is

-14-

within the range of about 34 to about 38 mmHg, and the  $pO_2$  is within the range of about 83 to about 89 mmHg. In a most preferred embodiment, following equilibration, the solution will have a pH of about 7.4, a  $pCO_2$  of about 36 mmHg and a  $pO_2$  of about 86 mmHg at 37°C.

In another preferred embodiment, the solution is representative of the blood of a human patient experiencing acidosis. In this embodiment, the solution has the following composition:

	<u>Component</u>	<u>Range</u>	<u>Specific</u>
	Total Bovine		
	Hemoglobin	8-10 g/dl	9 g/dl
	2-Phenoxyethanol	0.09-0.11 g/dl	0.1 g/dl
15	Benzethonium Chloride	0.09-0.11 g/dl	0.1 g/dl
	Sodium Selenite	1-20 uM	2 uM
	Phytic Acid	0.1-1.4 mM	0.14 mM
	Reduced Glutathione	0.5-8 mM	0.8 mM
	Sodium Bicarbonate	19.5-21.5 mM	20.5 mM
20	HEPES	20-40 mM	25 mM.

Upon use, this solution is equilibrated via tonometry with a gas mixture which comprises by volume 7%  $CO_2$ , 7%  $O_2$  and 86%  $N_2$ . This equilibration is carried out at 37°C until the pH is within the range of about 7.12 to about 7.18, the  $pCO_2$  is within the range of about 48 to about 52 mmHg, and the  $pO_2$  is within the range of about 47 to about 53 mmHg. In a most preferred embodiment, following equilibration, the solution will have a pH of about 7.15, a  $pCO_2$  of about 50 mmHg and a  $pO_2$  of about 50 mmHg at 37°C.



-15-

In yet another preferred embodiment, the solution is representative of the blood of a human patient experiencing alkalosis. In this embodiment, the solution has the following composition:

05	<u>Component</u>	<u>Range</u>	<u>Specific</u>
	Total Bovine		
	Hemoglobin	17-19 g/dl	18 g/dl
	2-Phenoxyethanol	0.09-0.11 g/dl	0.1 g/dl
	Benzethonium Chloride	0.09-0.11 g/dl	0.1 g/dl
10	Sodium Selenite	1-20 uM	2 uM
	Phytic Acid	0.1-1.4 mM	0.14 mM
	Reduced Glutathione	0.5-8 mM	0.8 mM
	Sodium Bicarbonate	19-21 mM	19.82 mM
	HEPPSO	20-40 mM	25 mM

15        Upon use, this solution is equilibrated via tonometry with a gas mixture which comprises by volume 2.8 CO<sub>2</sub>, 20% O<sub>2</sub> and 77.2% N<sub>2</sub>. This equilibration is carried out at 37°C until the pH is within the range of about 7.58 to about 7.62, the

20        pCO<sub>2</sub> is within the range of about 19 to about 21 mmHg, and the pO<sub>2</sub> is within the range of about 140 to about 144 mmHg. In a most preferred embodiment, the solution will have a pH of about 7.6, a pCO<sub>2</sub> of about 20 mmHg and a pO<sub>2</sub> of about 142 mmHg at 37°C.

25        A preferred batch protocol for producing a large volume of reference solution representative of normal blood is given below:

Reagents:

(A) Bovine red blood cell suspension

-16-

- (B) Aqueous solution containing 2% each of  
2-Phenoxyethanol and Benzethonium Chloride
- (C) Aqueous solution containing 0.072%  $\text{Na}_2\text{SeO}_3$   
and 1.87% Phytic Acid
- 05 (D) 119.15 g HEPES plus 11.03 g NaOH in 1  
liter of deionized water to provide a net  
pH of about 7.4
- (E) 4.92 g reduced glutathione in 1 liter  
deionized water

10 Procedures:

- 1) Wash the bovine red blood cell suspension  
with saline solution at least four times.
- 2) Adjust the total hemoglobin content to  
18.0 g% using deionized water, (some red  
cell lysis may occur at this point).
- 15 3) Complete the red cell lysis by adding 50  
ml of reagent (B) to 800 ml of the  
partially lysed red cell suspension,  
adding a small amount at a time and mixing  
thoroughly to avoid hemoglobin aggregation.
- 20 4) Add 50 ml each of reagents (C) and (D) to  
850 ml of the completely lysed hemoglobin  
solution obtained in step (3).
- 5) Add 50 ml of reagent (E) to bring the  
total volume to one liter.
- 25 6) Adjust the pH of the solution to 7.4 by  
slowly adding a 4% aqueous NaOH solution  
under a nitrogen (100%) atmosphere. (The  
total volume of NaOH added should be
- 30 recorded to determine whether addition of

-17-

reagents previously described is necessary.)

- 7) Add 2.45 g/l solid  $\text{NaHCO}_3$  to the solution.
- 05 8) Filter the hemoglobin solution through a filter having a maximum pore size of about 0.2  $\mu\text{m}$ .
- 9) Add a small quantity of a silicone-based defoamer such as antifoam A concentrate to the solution.
- 10 10) Flush the hemoglobin solution with nitrogen gas containing about 4 to 7%  $\text{CO}_2$  at room temperature until the oxyhemoglobin value is less than 1% of the total hemoglobin.
- 15 11) Subdivide the solution into individual sealed ampules having the gas mixture of step 10) contained in the head space.

The resulting, sealed ampules of reference solution can be stored at room temperature for a period of at least about three months and at 4°C for at least about one year. When used for monitoring of blood-gas/CO-Oximetry equipment is desired, the ampule is opened and subjected to tonometry with a predetermined gas mixture as previously described.

25 The resulting reference solution has many characteristics of fresh, human whole blood including physiological ranges of pH,  $\text{pCO}_2$ ,  $\text{pO}_2$ , total hemoglobin content and fractions of hemoglobin species.

-18-

Equivalents

Those skilled in the art will recognize, or be able to ascertain using no more than routine experimentation, many equivalents to the specific embodiments of the invention described herein. Such equivalents are intended to be covered by the following claims.

-19-

CLAIMS

- 05 1. A stable hemoglobin reference solution comprising an aqueous solution containing:
- a) hemoglobin;
  - b) a source of bicarbonate ions;
  - c) a source of thio functional groups;
  - d) a source of metal ion catalyst;
  - e) at least one organic polyphosphate; and,
  - f) an organic buffer.
- 
- 10 2. A solution of Claim 1 wherein the hemoglobin comprises 95% reduced hemoglobin.
- 15 3. A solution of Claim 2 wherein the solution contains:
- a) less than about 3%, by weight, methemoglobin,
  - b) less than about 2%, by weight, carboxyhemoglobin,
  - c) less than about 2%, by weight, oxyhemoglobin, and
  - 20 d) less than about 0.3%, by volume, oxygen.
4. A stable reference solution of Claim 3 additionally including an antimicrobial agent.
- 25 5. A stable reference solution of in Claim 4 having a  $P_{50}$  value of between about 18 to about 34 mmHg after tonometry with an appropriate gas mixture of  $O_2$ ,  $CO_2$  and  $N_2$ .

-20-

6. A stable reference solution of Claim 5 wherein the hemoglobin comprises mammalian hemoglobin.
7. A stable reference solution of Claim 6 wherein the hemoglobin comprises bovine hemoglobin.
- 05 8. A stable reference solution of Claim 7 wherein thio functional groups are provided by a sulfhydro compound selected from reduced glutathione, mercaptoethanol, cysteine or derivatives thereof.
- 10 9. A stable reference solution of Claim 8 wherein the source of metal ion catalyst comprises a salt of metals selected from the group consisting of selenium, tellurium and copper.
- 15 10. A stable reference solution of Claim 9 wherein the organic polyphosphate compound is selected from the group consisting of 2,3-diphosphoglycerate, phytic acid and adenosine triphosphate.
- 20 11. A stable reference solution of Claim 1 contained in a sealed container having a head space containing CO<sub>2</sub> and a gas which is inert to the reference solution.

-21-

12. A stable reference solution which comprises an aqueous solution containing:
- 05 a) between about 13 and about 15 g/dl bovine hemoglobin solution, said hemoglobin solution comprising fractions of:
- 10 i) at least about 95% reduced hemo-  
globin;  
ii) less than about 3% methemoglobin;  
iii) less than about 2% carboxyhemoglobin;  
15 iv) less than about 2% oxyhemoglobin; and  
v) less than about 0.3% by volume  
oxygen;
- 20 b) between about 0.09 and about 0.11 g/dl 2-phenoxyethanol;
- 15 c) between about 0.09 and about 0.11 g/dl benzethonium chloride;
- 20 d) a concentration of sodium selenite ranging from about 1.0 to about 20  $\mu$ M;
- 25 e) a concentration of phytic acid ranging from about 0.1 to about 1.4 mM;
- 20 f) a concentration of reduced glutathione ranging from about 0.5 to about 8 mM;
- 25 g) a concentration of sodium bicarbonate ranging from about 28.5 to about 30.5 mM;  
and
- h) a concentration of HEPES buffer ranging from about 20 to about 40 mM.
13. A stable reference of Claim 12 contained within a sealed container, said container having a head space containing at least about 93% inert gas selected from the group consisting of nitrogen, noble gasses and mixtures thereof, and the remainder of the gas comprising CO<sub>2</sub>.
- 30

-22-

14. A stable reference solution which comprises an aqueous solution containing:
- 05 a) between about 8 and about 10 g/dl bovine hemoglobin solution, said hemoglobin solution comprising fractions of:
- 10 i) at least about 95% reduced hemoglobin,
- ii) less than about 3% methemoglobin,
- iii) less than about 2% carboxyhemoglobin,
- iv) less than about 2% oxyhemoglobin, and
- 15 v) less than about 0.3% by volume oxygen;
- b) between about 0.09 and about 0.11 g/dl 2-phenoxyethanol;
- 20 c) between about 0.09 and about 0.11 g/dl benzethonium chloride;
- d) a concentration of sodium selenite ranging from about 1 to about 20  $\mu$ M;
- e) a concentration of phytic acid ranging from about 0.1 to about 1.4 mM;
- 25 f) a concentration of reduced glutathione ranging from about 0.5 to about 8 mM;
- g) a concentration of sodium bicarbonate ranging from about 19.5 to about 21.5 mM; and
- h) a concentration of HEPES buffer ranging from about 20 to about 40 mM.



-23-

15. A stable reference solution of Claim 14 contained within a sealed container, said container having a head space containing at least about 93% inert gas selected from the group consisting of nitrogen, noble gasses and mixtures thereof, and the remainder of the gas comprising CO<sub>2</sub>.
- 05
16. A stable reference solution for which comprises an aqueous mixture containing:
- 10 a) between about 17 and about 19 g/dl bovine hemoglobin solution, said hemoglobin solution comprising fractions of:
- 15 i) at least about 95% reduced hemoglobin;
- ii) less than about 3% methemoglobin;
- iii) less than about 2% carboxyhemoglobin;
- iv) less than about 2% oxyhemoglobin;
- v) less than about 0.3% by volume oxygen;
- 20 b) between about 0.09 and about 0.11 g/dl 2-phenoxyethanol;
- c) between about 0.09 and about 0.11 g/dl benzethonium chloride;
- 25 d) a concentration of sodium selenite ranging from about 1 to about 20  $\mu$ M;
- e) a concentration of phytic acid ranging from about 0.1 to about 1.4 mM;
- f) a concentration of reduced glutathione ranging from about 0.5 to about 8 mM;


-24-

- g) a concentration of sodium bicarbonate ranging from about 19 to about 21 mM; and
- h) a concentration of HEPPSO buffer ranging from about 20 to about 40 mM.

- 05 17. A stable reference solution of Claim 16 contained within a sealed container, said container having a head space containing at least about 93% inert gas selected from the group consisting of nitrogen, noble gasses and mixtures thereof, and the remainder of the gas comprising CO<sub>2</sub>.
- 10

# INTERNATIONAL SEARCH REPORT

International Application No PCT/US 90/01422

<b>I. CLASSIFICATION OF SUBJECT MATTER</b> (if several classification symbols apply, indicate all) <sup>6</sup>		
According to International Patent Classification (IPC) or to both National Classification and IPC		
IPC5: G 01 N 33/72, G 01 N 33/96		
<b>II. FIELDS SEARCHED</b>		
Minimum Documentation Searched <sup>7</sup>		
Classification System	Classification Symbols	
IPC5	G 01 N	
Documentation Searched other than Minimum Documentation to the Extent that such Documents are Included in Fields Searched <sup>8</sup>		
<b>III. DOCUMENTS CONSIDERED TO BE RELEVANT<sup>9</sup></b>		
Category *	Citation of Document, <sup>11</sup> with indication, where appropriate, of the relevant passages <sup>12</sup>	Relevant to Claim No. <sup>13</sup>
X	EP, A1, 0185404 (STICHTING GASTRANSPORT) 25 June 1986, see the whole document  --	1-17
X	EP, A1, 0181033 (AKZO N.V.) 14 May 1986, see the whole document  --	1-17
A	US, A, 4485174 (CHIANG ET AL) 27 November 1984, see the whole document  --	1-17
A	US, A, 4469792 (SIMMONDS ET AL) 4 September 1984, see the whole document  --  -----	1-17
<p><sup>a</sup> Special categories of cited documents: <sup>10</sup></p> <p>"A" document defining the general state of the art which is not considered to be of particular relevance</p> <p>"E" earlier document but published on or after the international filing date</p> <p>"L" document which may throw doubts on priority claim(s) or which is cited to establish the publication date of another citation or other special reason (as specified)</p> <p>"O" document referring to an oral disclosure, use, exhibition or other means</p> <p>"P" document published prior to the international filing date but later than the priority date claimed</p> <p>"T" later document published after the international filing date or priority date and not in conflict with the application but cited to understand the principle or theory underlying the invention</p> <p>"X" document of particular relevance, the claimed invention cannot be considered novel or cannot be considered to involve an inventive step</p> <p>"Y" document of particular relevance, the claimed invention cannot be considered to involve an inventive step when the document is combined with one or more other such documents, such combination being obvious to a person skilled in the art.</p> <p>"&amp;" document member of the same patent family</p>		
<b>IV. CERTIFICATION</b>		
Date of the Actual Completion of the International Search	Date of Mailing of this International Search Report	
13th June 1990	28.06.90	
International Searching Authority	Signature of Authorized Officer	
EUROPEAN PATENT OFFICE	 T. TAZELAAR	

**ANNEX TO THE INTERNATIONAL SEARCH REPORT  
ON INTERNATIONAL PATENT APPLICATION NO. PCT/US 90/01422**

SA 35545

This annex lists the patent family members relating to the patent documents cited in the above-mentioned international search report.  
The members are as contained in the European Patent Office EDP file on 07/05/90  
The European Patent office is in no way liable for these particulars which are merely given for the purpose of information.

Patent document cited in search report	Publication date	Patent family member(s)	Publication date
EP-A1- 0185404	25/06/86	NL-A- 8403425	02/06/86
EP-A1- 0181033	14/05/86	AU-B- 589582	19/10/89
		AU-D- 4923785	15/05/86
		CA-A- 1246429	13/12/88
		JP-A- 61114168	31/05/86
		US-A- 4711852	08/12/87
US-A- 4485174	27/11/84	EP-A-B- 0158684	23/10/85
		JP-A- 60231166	16/11/85
US-A- 4469792	04/09/84	US-A- 4401652	30/08/83

For more details about this annex : see Official Journal of the European patent Office, No. 12/82